

## XP-002256825

AN - 1997-095485 [09]  
AP - JP19950144406 19950612  
CPY - MITV  
DC - B04  
FS - CPI  
IC - C07B63/00 ; C07K1/18 ; C07K1/22  
MC - B04-H19 B04-L01 B04-N04  
M1 - [01] M423 M720 M903 N164 V752  
M2 - [02] A220 A940 B115 B701 B713 B720 B815 B831 C101 C108 C550 C730 C802  
C804 C805 C807 M411 M781 M903 M904; R03521-U  
PA - (MITV ) MITSUBISHI MATERIALS CORP  
PN - JP8333387 A 19961217 DW199709 C07K1/18 004pp  
PR - JP19950144406 19950612  
XA - C1997-030580  
XIC - C07B-063/00 ; C07K-001/18 ; C07K-001/22  
AB - J08333387 Isolation and purificn. of protein by ion exchange  
absorption and sepn., in which a protein contg. soln. is contacted  
with hydroxyapatite treated with a metal salt or phosphate soln., to  
selectively absorb and separate acidic or basic protein.  
- The hydroxyapatite prepd. by baking at temps. of 1,000 deg.C or lower.  
Hydroxyapatite is baked at temps. of at most 1,000 deg.C, pref.  
800-1,000 deg.C and pulverised to 5-500 micron, particles. Then, the  
baked hydroxyapatite is contacted with a metal salt soln. (e.g. Ca,  
Mg, Na, K, Mn, Li, Ni, Cu, Zn, Fe and Ag, pref. Ca and Mg ions) at  
concns. of 2 mM or less, pref. 1-2 mM.  
- ADVANTAGE - Selective and ion exchange absorption and sepn. of acidic  
or basic protein by simple method.  
- In an example, hydroxyapatite particles having 50-150 micron sizes  
were baked at 900 deg.C and filled in a column of 0.9 cm diameter and  
15 cm length. The column was equilibrated with 30 mM potassium  
phosphate buffer and a mixt. of albumin (BSA), fibrinogen, lysozyme  
and cytochrome having isoelectric point (pI) of 4.7, 5.6, 10.8 and  
10.6, respectively, was poured. Then, the column was linear gradiently  
developed with 0.03-0.35M potassium phosphate buffer to selectively  
absorb lysozyme and cytochrome without absorbing albumin or fibrinogen.  
- (Dwg.1/2)  
CN - R03521-U  
IW - ISOLATE PURIFICATION PROTEIN ION EXCHANGE ABSORB HYDROXY APATITE  
IKW - ISOLATE PURIFICATION PROTEIN ION EXCHANGE ABSORB HYDROXY APATITE  
NC - 001  
OPD - 1995-06-12  
ORD - 1996-12-17  
PAW - (MITV ) MITSUBISHI MATERIALS CORP  
TI - Isolation and purificn. of protein - by ion exchange and absorption  
using hydroxy:apatite